


Anti-Ep400 Polyclonal Antibody For Wb, Elisa - Q96L91

Item Number: CM0014220



Introduction

Rabbit polyclonal antibody targeting human EP400 (p400), a 343 kDa NuA4 complex component. Validated for WB and ELISA applications. Ideal for chromatin remodeling and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	EP400 Rabbit pAb
Remarks/Alias	P400; CAGH32; TNRC12; EP400
Species	Human
GeneID (Human)	57634
GeneID	57634
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96L91
Protein Weight	343kDa
Shipping	Ice bag