

Anti-Epac1/Rapgef3 Polyclonal Antibody For Wb, Elisa - 095398

Item Number: CM0019966



Introduction

Rabbit polyclonal antibody targeting human Epac1/RAPGEF3, a cAMP-regulated Rap guanine nucleotide exchange factor. Validated for WB and ELISA; cross-reacts with human and mouse epitopes. Predicted molecular weight: 104 kDa.

[Learn More](#)

Parameter	Value
Product Name	Epac1/RAPGEF3 Rabbit pAb
Remarks/Alias	EPAC; EPAC1; bcm910; HSU79275; CAMP-GEFI; Epac1/RAPGEF3
Species	Human
Gene ID (Human)	10411
Gene ID	10411
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Epac1/Epac1/RAPGEF3 (NP_001092002.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95398
Protein Weight	104kDa
Shipping	Ice bag