

Anti-Epac1/Rapgef3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O95398

Item Number: CM0009390



Introduction

Rabbit monoclonal antibody against human Epac1/RAPGEF3 (SWISS: O95398) for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Epac1 is a cAMP-activated GEF for Rap1, involved in angiogenesis and barrier function.

[Learn More](#)

Parameter	Value
Product Name	Epac1/RAPGEF3 Rabbit mAb
Remarks/Alias	EPAC; Epac1/RAPGEF3
Species	Human
Gene ID (Human)	10411
Gene ID	10411
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 824-923 of human Epac1/RAPGEF3 (O95398).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95398
Protein Weight	104kDa
Shipping	Ice bag