

Anti-Epas1/Hif2A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99814

Item Number: CM0030475



Introduction

High-quality rabbit polyclonal antibody against human EPAS1/HIF2 α (SWISS: Q99814). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for hypoxia and angiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	EPAS1/HIF2 α Rabbit pAb
Remarks/Alias	HLF; MOP2; ECT4; HIF2A; PASD2; bHLHe73; EPAS1/HIF2 α
Species	Human
GeneID (Human)	2034
GeneID	2034
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99814
Protein Weight	96kDa
Shipping	Ice bag