


Anti-Epb41 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11171

Item Number: CM0005114



Introduction

Recombinant rabbit monoclonal antibody targeting human EPB41 (4.1R), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for membrane cytoskeleton, erythrocyte biology, and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	EPB41 Rabbit mAb
Remarks/Alias	HE; EL1; 4.1R
Species	Human
Gene ID (Human)	2035
Gene ID	2035
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 711-810 of human EPB41 (NP_001159477.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P11171
Protein Weight	97kDa
Shipping	Ice bag