

Anti-Epb41L2 Polyclonal Antibody For Wb, If/Icc, Elisa - O43491

Item Number: CM0029530

Introduction

Rabbit polyclonal antibody targeting human EPB41L2 (4.1G), validated for WB, IF/ICC, and ELISA. Suitable for studying cytoskeleton, mitosis, and cell cortex function.

[Learn More](#)

Parameter	Value
Product Name	EPB41L2 Rabbit pAb
Remarks/Alias	4.1G; 4.1-G; EPB41L2
Species	Human
GeneID (Human)	2037
GeneID	2037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human EPB41L2 (NP_001186317.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O43491
Protein Weight	113kDa
Shipping	Ice bag