

Anti-Epb42 (Erythrocyte Protein 4.2) Rabbit Polyclonal Antibody For If/Icc And Elisa - P16452

Item Number: CM0028927



Introduction

EPB42 Rabbit pAb is a polyclonal antibody targeting erythrocyte protein 4.2, validated for IF/ICC and ELISA in human and mouse. Ideal for studying erythrocyte membrane stability, anemia, and related disorders.

[Learn More](#)

Parameter	Value
Product Name	EPB42 Rabbit pAb
Remarks/Alias	PA; SPH5
Species	Human
Gene ID (Human)	2038
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 589-691 of human EPB42 (NP_001107606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16452
Protein Weight	77kDa
Shipping	Ice bag