

Anti-Epcam Polyclonal Antibody For Wb, Elisa - Q99Jw5

Item Number: CM0023064

Introduction

Rabbit polyclonal anti-Epcam antibody for WB and ELISA. Reacts with human, mouse, rat. Targets epithelial cell adhesion molecule (CD326, EGP314). Immunogen: recombinant mouse Epcam (aa 24-266).

[Learn More](#)

Parameter	Value
Product Name	Epcam Rabbit pAb
Remarks/Alias	EGP; Ly74; gp40; CD326; EGP-2; TROP1; Egp314; Ep-CAM; EpCAM1; Tacsd1; GA733-2; Tacstd1; Epcam
Species	Mouse
Gene ID	17075
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-266 of mouse Epcam (NP_032558.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99JW5
Protein Weight	35kDa
Shipping	Ice bag