

Anti-Epcam Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P16422

Item Number: CM0002575



Introduction

Rabbit monoclonal antibody against human EpCAM (CD326) for Western blot, IHC-P, and ELISA. Recognizes human, mouse, and rat EpCAM. Useful for epithelial cell adhesion research and cancer stem cell studies.

[Learn More](#)

Parameter	Value
Product Name	EpCAM Rabbit mAb
Remarks/Alias	ESA; KSA; M4S1; MK-1; DIAR5; EGP-2; EGP40; KS1/4; MIC18; TROP1; BerEp4; EGP314; HNPCC8; LYNCH8; MOC-31; Ber-Ep4; TACSTD1; EpCAM
Species	Human
Gene ID (Human)	4072
Gene ID	4072
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human EpCAM (P16422))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16422
Protein Weight	35kDa
Shipping	Ice bag