

## Anti-Epcam Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P16422

Item Number: CM0014628



### Introduction

EpCAM rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for epithelial cell adhesion molecule research.

[Learn More](#)

| Parameter        | Value                                                                                                                           |
|------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EpCAM Rabbit pAb                                                                                                                |
| Remarks/Alias    | ESA; KSA; M4S1; MK-1; DIAR5; EGP-2; EGP40; KS1/4; MIC18; TROP1; BerEp4; EGP314; HNPCC8; LYNCH8; MOC-31; Ber-Ep4; TACSTD1; EpCAM |
| Species          | Human                                                                                                                           |
| GeneID (Human)   | 4072                                                                                                                            |
| GeneID           | 4072                                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of human EpCAM (NP_002345.2)               |
| Source           | Rabbit                                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                               |
| SWISS            | P16422                                                                                                                          |
| Protein Weight   | 35kDa                                                                                                                           |
| Shipping         | Ice bag                                                                                                                         |