

Anti-Epcr/Cd201 Monoclonal Antibody For Wb, Elisa - Q64695

Item Number: CM0007318

Introduction

Recombinant rabbit monoclonal antibody targeting mouse EPCR/CD201. Validated for WB and ELISA. Involved in blood coagulation through protein C activation. Suitable for hemostasis and endothelial cell research.

[Learn More](#)

Parameter	Value
Product Name	EPCR/CD201 Rabbit mAb
Remarks/Alias	Ccca; Epcr; Ccd41; EPCR/CD201
Species	Mouse
Gene ID (Human)	10544
Gene ID	19124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of mouse EPCR/CD201 (NP_035301.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q64695
Protein Weight	27kDa
Shipping	Ice bag