


Anti-Epcr/Cd201 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q9Unn8

Item Number: CM0004436



Introduction

Rabbit monoclonal antibody targeting human EPCR/CD201 protein, validated for WB, FC, and ELISA applications, with cross-reactivity to mouse. Suitable for coagulation and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	EPCR/CD201 Rabbit mAb
Remarks/Alias	PROCR; CCCA; CCD41; EPCR; protein C receptor; EPCR/CD201
Species	Human
GeneID (Human)	10544
GeneID	10544
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-210 of human CD201/EPCR (NP_006395.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNN8
Protein Weight	26kDa
Shipping	Ice bag