


Anti-Eph Receptor B3 (Ephb3) Monoclonal Antibody For Wb, Elisa - P54753

Item Number: CM0003021



Introduction

Rabbit monoclonal antibody against human Eph receptor B3 (EPHB3), validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Targets a key receptor tyrosine kinase involved in axon guidance and synapse development.

[Learn More](#)

Parameter	Value
Product Name	Eph receptor B3 (EPHB3) Rabbit mAb
Remarks/Alias	EK2; ETK2; HEK2; TYRO6; Eph receptor B3 (EPHB3)
Species	Human
GeneID	2049 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human Eph receptor B3 (EPHB3) (P54753)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54753
Protein Weight	110kDa
Shipping	Ice bag