

Anti-Epha2 Rabbit Monoclonal Antibody For Wb And Elisa - P29317

Item Number: CM0005039



Introduction

Rabbit monoclonal antibody targeting human EPHA2 (Ephrin type-A receptor 2). Validated for Western blot and ELISA. Suitable for research on cell adhesion, migration, differentiation, and viral entry.

[Learn More](#)

Parameter	Value
Product Name	EPHA2 Rabbit mAb
Remarks/Alias	ECK; CTPA; ARCC2; CTPP1; CTRCT6; EPHA2
Species	Human
GeneID (Human)	1969
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-534 of human EPHA2. (NP_004422.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P29317
Protein Weight	108kDa
Shipping	Ice bag