

Anti-Epha3 Rabbit Polyclonal Antibody For Wb, Elisa - P29320

Item Number: CM0027985



Introduction

Anti-EPHA3 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat EPHA3 (110 kDa). Suitable for studies of ephrin receptor signaling, cell adhesion, and developmental processes.

[Learn More](#)

Parameter	Value
Product Name	EPHA3 Rabbit pAb
Remarks/Alias	EK4; ETK; HEK; ETK1; HEK4; TYRO4; EPHA3
Species	Human
GeneID (Human)	2042
GeneID	2042
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 400-545 of human EPHA3 (NP_005224.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29320
Protein Weight	110kDa
Shipping	Ice bag