

Anti-Epha3/Epha4/Epha5 Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P29320 / P54764 / P54756

Item Number: CM0019401



Introduction

Rabbit polyclonal antibody targeting human EPHA3, EPHA4, and EPHA5 receptors. Validated for IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Suitable for studying ephrin receptor signaling and cell adhesion in development and disease research.

[Learn More](#)

Parameter	Value
Product Name	EPHA3/EPHA4/EPHA5 pAb
Species	Human
Gene ID (Human)	2042/2043/2044
Gene ID	2042/2043/2044
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-800 of human EPHA3/EPHA4/EPHA5 (NP_005224.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29320/P54764/P54756
Shipping	Ice bag