

Anti-Epha4 Polyclonal Antibody For Wb, If/Icc, Elisa - P54764

Item Number: CM0029427



Introduction

Rabbit polyclonal antibody against EPHA4 (EK8/SEK/TYRO1) validated for Human, Mouse, Rat in WB, IF/ICC, ELISA. Targets ephrin receptor involved in nervous system development and repair.

[Learn More](#)

Parameter	Value
Product Name	EPHA4 Rabbit pAb
Remarks/Alias	EK8; SEK; HEK8; TYRO1; EPHA4
Species	Human
GeneID (Human)	2043
GeneID	2043
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-100 of human EPHA4 (NP_004429.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54764
Protein Weight	110kDa
Shipping	Ice bag