

Anti-Ephb2 Polyclonal Antibody For Wb, Elisa - P29323

Item Number: CM0024712



Introduction

Rabbit polyclonal antibody targeting human EPHB2 (NP_059145.2, amino acids 250-540). Validated for WB and ELISA; cross-reacts with human and mouse. Ideal for developmental neuroscience, angiogenesis, and tumor suppression research.

[Learn More](#)

Parameter	Value
Product Name	EPHB2 Rabbit pAb
Remarks/Alias	DRT; EK5; ERK; CAPB; Hek5; PCBC; EPHT3; Tyro5; BDPLT22; EPHB2
Species	Human
GeneID (Human)	2048
GeneID	2048
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 250-540 of human EPHB2 (NP_059145.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29323
Protein Weight	117kDa
Shipping	Ice bag