


Anti-Ephb4 Rabbit Polyclonal Antibody For Wb, Elisa - P54760

Item Number: CM0029458



Introduction

Rabbit polyclonal antibody against human EPHB4, suitable for WB and ELISA, with cross-reactivity to mouse. Designed for research on angiogenesis, cardiovascular development, and cell migration using recombinant protein immunogen.

[Learn More](#)

Parameter	Value
Product Name	EPHB4 Rabbit pAb
Remarks/Alias	HTK; MYK1; HFASD; CMAVM2; LMPHM7; TYRO11; EPHB4
Species	Human
Gene ID (Human)	2050
Gene ID	2050
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54760
Protein Weight	108kDa
Shipping	Ice bag