

Anti-Ephb6 Polyclonal Antibody For Wb, Elisa - O15197

Item Number: CM0017536



Introduction

Anti-EPHB6 polyclonal antibody (rabbit) for WB, ELISA. Reacts with human, mouse. Immunogen: recombinant fusion protein covering amino acids 615-729. Target: kinase-defective ephrin receptor involved in cell adhesion and migration. SWISS: O15197.

[Learn More](#)

Parameter	Value
Product Name	EPHB6 Rabbit pAb
Remarks/Alias	HEP; EPHB6
Species	Human
GeneID (Human)	2051
GeneID	2051
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 615-729 of human EPHB6 (NP_001267724.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15197
Protein Weight	111kDa
Shipping	Ice bag