

Anti-Ephrin B2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P52799

Item Number: CM0003489



Introduction

High-quality rabbit monoclonal antibody against human Ephrin B2 (EFNB2). Validated for WB, IHC-P, ELISA with cross-reactivity in human, mouse, rat. Ideal for angiogenesis, neurodevelopment, and viral entry studies.

[Learn More](#)

Parameter	Value
Product Name	Ephrin B2 Rabbit mAb
Remarks/Alias	HTKL; EPLG5; Htk-L; LERK5; ephrin-B2; Ephrin B2
Species	Human
Gene ID (Human)	1948
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human EFNB2 (P52799)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52799
Protein Weight	37 kDa
Shipping	Ice bag