

Anti-Ephrin B2 Rabbit Polyclonal Antibody For Wb, Elisa - P52799

Item Number: CM0015253



Introduction

High-quality rabbit polyclonal antibody against human Ephrin B2, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying neuronal development and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	Ephrin B2 Rabbit pAb
Remarks/Alias	HTKL; EPLG5; Htk-L; LERK5; ephrin-B2; Ephrin B2
Species	Human
GeneID (human)	1948
GeneID	1948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-229 of human Ephrin B2 (NP_004084.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P52799
Protein Weight	37kDa
Shipping	Ice bag