

Anti-Ephx1 Rabbit Polyclonal Antibody For Wb, Elisa - P07099

Item Number: CM0023744

Introduction

High-quality rabbit polyclonal antibody targeting EPHX1 (epoxide hydrolase 1, mEH). Validated for WB and ELISA applications, cross-reacts with human and mouse. Ideal for metabolism and detoxification studies.

[Learn More](#)

Parameter	Value
Product Name	EPHX1 Rabbit pAb
Remarks/Alias	MEH; EPHX; EPOX; HYL1; EPHX1
Species	Human
GeneID	2052
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human EPHX1 (NP_000111.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07099
Protein Weight	53kDa
Shipping	Ice bag