


Anti-Ephx2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P34913

Item Number: CM0021756



Introduction

Rabbit polyclonal EPHX2 antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat EPHX2. Ideal for investigating epoxide hydrolase and lipid phosphatase functions in detoxification and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	EPHX2 Rabbit pAb
Remarks/Alias	CEH; SEH; ABHD20; EPHX2
Species	Human
GeneID (Human)	2053
GeneID	2053
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-240 of human EPHX2 (NP_001970.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P34913
Protein Weight	63kDa
Shipping	Ice bag