


Anti-Epm2A Rabbit Polyclonal Antibody For Wb, Elisa - 095278

Item Number: CM0028396



Introduction

Rabbit polyclonal antibody targeting human EPM2A/laforin, validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for Lafora disease, glycogen metabolism, and neurology research.

[Learn More](#)

Parameter	Value
Product Name	EPM2A Rabbit pAb
Remarks/Alias	EPM2; MELF; EPM2A
Species	Human
Gene ID (Human)	7957
Gene ID	7957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 244-331 of human EPM2A (NP_005661.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95278
Protein Weight	37kDa
Shipping	Ice bag