


Anti-Epn1 (Epsin-1) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6I3

Item Number: CM0024065



Introduction

Validated rabbit polyclonal antibody against human EPN1 (Epsin-1). Suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Target involved in clathrin-mediated endocytosis and membrane curvature.

[Learn More](#)

Parameter	Value
Product Name	EPN1 Rabbit pAb
Remarks/Alias	EPN1
Species	Human
GeneID (human)	29924
GeneID	29924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-180 of human EPN1 (NP_037465.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6I3
Protein Weight	60kDa
Shipping	Ice bag