


## Anti-Epo Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01588

Item Number: CM0027795



### Introduction

Rabbit polyclonal antibody targeting human Erythropoietin (EPO). Validated for WB, IF-P, IHC-P, ELISA, and cross-reacts with mouse and rat. Ideal for studying erythropoiesis and related disorders.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EPO Rabbit pAb                                                                                                         |
| Remarks/Alias    | EP; DBAL; ECTY5; MVCD2; EPO                                                                                            |
| Species          | Human                                                                                                                  |
| Gene ID          | 2056                                                                                                                   |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 93-193 of human EPO (NP_000790.2) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, IF-P, IHC-P, ELISA                                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                                      |
| SWISS            | P01588                                                                                                                 |
| Protein Weight   | 21kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |