


## Anti-Epor Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P19235

Item Number: CM0028876



### Introduction

Rabbit polyclonal antibody targeting human erythropoietin receptor (EPOR), validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for erythropoietin signaling and erythroid differentiation studies.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | EPOR Rabbit pAb                                                                                         |
| Remarks/Alias    | EPO-R; EPOR                                                                                             |
| Species          | Human                                                                                                   |
| Gene ID (Human)  | 2057                                                                                                    |
| Gene ID          | 2057                                                                                                    |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 274-508 of human EPOR (NP_000112.1). |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                 |
| Cross Reactivity | Human, Mouse, Rat                                                                                       |
| SWISS            | P19235                                                                                                  |
| Protein Weight   | 55kDa                                                                                                   |
| Shipping         | Ice bag                                                                                                 |