

Anti-Eppin Polyclonal Antibody For Wb, Elisa - 095925

Item Number: CM0023387

Introduction

Anti-EPPIN rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat. EPPIN is a serine protease inhibitor critical for male fertility and sperm motility. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	EPPIN Rabbit pAb
Remarks/Alias	CT71, CT72, WAP7, WFDC7, SPINLW1, dJ461P17.2, EPPIN
Species	Human
GeneID (human)	57119
GeneID	57119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of human EPPIN (NP_065131.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95925
Protein Weight	15kDa
Shipping	Ice bag