

Anti-Eprs Rabbit Polyclonal Antibody For Wb, Elisa - P07814

Item Number: CM0013766



Introduction

Rabbit polyclonal antibody against human EPRS (bifunctional glutamate/proline--tRNA ligase), validated for WB and ELISA. Suitable for research in translation regulation, aminoacyl-tRNA synthetase complex, and GAIT complex studies.

[Learn More](#)

Parameter	Value
Product Name	EPRS Rabbit pAb
Remarks/Alias	EARS; EPRS; PARS; QARS; QPRS; HLD15; PIG32; GLUPRORS
Species	Human
Gene ID (Human)	2058
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07814
Protein Weight	171kDa
Shipping	Ice bag