


Anti-Eps15 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P42566

Item Number: CM0024967



Introduction

Rabbit polyclonal antibody targeting human EPS15 (EGFR substrate 15, AF-1p), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat. Suitable for endocytosis and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	EPS15 Rabbit Polyclonal Antibody
Remarks/Alias	AF1P; AF-1P; MLLT5; EPS15
Species	Human
Gene ID	2060
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 660-840 of human EPS15 (NP_001972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42566
Protein Weight	99kDa
Shipping	Ice bag