

Anti-Eps15L1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ubc2

Item Number: CM0029167

Introduction

CamelBio's anti-EPS15L1 rabbit polyclonal antibody targets human EPS15L1 (Eps15R). Validated for Western blot and ELISA, it supports research into receptor-mediated endocytosis and integrin trafficking.

[Learn More](#)

Parameter	Value
Product Name	EPS15L1 Rabbit pAb
Remarks/Alias	EPS15R; EPS15L1
Species	Human
Gene ID (Human)	58513
Gene ID	58513
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 475-754 of human EPS15L1 (NP_001245304.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBC2
Protein Weight	94kDa
Shipping	Ice bag