


Anti-Eps8 Polyclonal Antibody For Wb, Elisa - Q12929

Item Number: CM0015970



Introduction

Rabbit polyclonal antibody targeting human EPS8, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying actin cytoskeleton dynamics, cell migration, and cancer.

[Learn More](#)

Parameter	Value
Product Name	EPS8 Rabbit pAb
Remarks/Alias	DFNB102; EPS8
Species	Human
GeneID (human)	2059
GeneID	2059
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human EPS8 (NP_004438.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12929
Protein Weight	92kDa
Shipping	Ice bag