


Anti-Erab/Hsd17B10 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99714

Item Number: CM0028124



Introduction

High-quality rabbit polyclonal antibody against human ERAB/HSD17B10, suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial metabolism and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	ERAB/HSD17B10 Rabbit pAb
Remarks/Alias	ABAD; CAMR; ERAB; HCD2; MHBD; HADH2; MRPP2; MRX17; MRX31; SCHAD; MRXS10; SDR5C1; HSD10MD; 17b-HSD10; DUPXp11.22; ERAB/HSD17B10
Species	Human
Gene ID (Human)	3028
Gene ID	3028
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99714
Protein Weight	27kDa
Shipping	Ice bag