


Anti-ErbB2Ip Rabbit Polyclonal Antibody For Wb, Elisa - Q96Rt1

Item Number: CM0028696



Introduction

High-quality rabbit polyclonal antibody targeting human ERBB2IP (Erbin, LAP2). Validated for WB and ELISA. Essential for studying ERBB2 signaling, keratinocyte differentiation, and NF-κB regulation in epithelia.

[Learn More](#)

Parameter	Value
Product Name	ERBB2IP Rabbit pAb
Remarks/Alias	LAP2; ERBB2IP; HEL-S-78
Species	Human
GeneID	55914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1152-1371 of human ERBB2IP (NP_061165.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96RT1
Protein Weight	158kDa
Shipping	Ice bag