


Anti-ErbB3/Her3 Monoclonal Antibody For Flow Cytometry - P21860

Item Number: CM0005349



Introduction

CamelBio offers -conjugated rabbit anti-human ErbB3/HER3 monoclonal antibody validated for flow cytometry. Targets the extracellular domain, ideal for studying HER3 signaling in human epithelial and brain tissues.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ErbB3/HER3 mAb
Remarks/Alias	ERBB3; ErbB-3; HER3; LCCS2; MDA-BF-1; c-erbB-3; c-erbB3; erbB3-S; p180-ErbB3; p45-sErbB3; p85-sErbB3; erb-b2 receptor tyrosine kinase 3
Species	Human
GeneID (Human)	2065
GeneID	2065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-643 of human ErbB3/HER3 (NP_001973.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21860
Protein Weight	148kDa
Shipping	Ice bag