


Anti-ErbB3/Her3 Monoclonal Antibody (, Rabbit) For Flow Cytometry - P21860

Item Number: CM0007754



Introduction

-conjugated rabbit monoclonal antibody targeting human ErbB3/HER3, validated for flow cytometry (FC). Suitable for detection of receptor tyrosine-protein kinase erbB-3 in cell surface and intracellular studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ErbB3/HER3 mAb
Remarks/Alias	HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3
Species	Human
GeneID (human)	2065
GeneID	2065
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-643 of human ErbB3/HER3 (NP_001973.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21860
Protein Weight	148kDa
Shipping	Ice bag