


Anti-ErbB3/Her3 Polyclonal Antibody For Wb, If/Icc, Elisa - P21860

Item Number: CM0017749



Introduction

Anti-ErbB3/HER3 polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat ErbB3/HER3. Targets a neuregulin receptor tyrosine kinase important in cancer research.

[Learn More](#)

Parameter	Value
Product Name	ErbB3/HER3 Rabbit pAb
Remarks/Alias	HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; ErbB3/HER3
Species	Human
Gene ID (Human)	2065
Gene ID	2065
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-285 of human ErbB3/HER3 (NP_001973.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21860
Protein Weight	148kDa
Shipping	Ice bag