


Anti-ErbB4/Her4 Monoclonal Antibody For Wb, Elisa - Q15303

Item Number: CM0002270



Introduction

Rabbit monoclonal antibody against human ErbB4/HER4, cross-reactive with mouse and rat. Suitable for WB and ELISA, targeting a receptor tyrosine kinase involved in development and cancer.

[Learn More](#)

Parameter	Value
Product Name	ErbB4/HER4 Rabbit mAb
Remarks/Alias	HER4; ALS19; p180erbB4; ErbB4/HER4
Species	Human
Gene ID (Human)	2066
Gene ID	2066
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1200-1308 of human ErbB4/HER4 (Q15303).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q15303
Protein Weight	147kDa
Shipping	Ice bag