


## Anti-ErbB4/Her4 Rabbit Polyclonal Antibody For Wb And Elisa - Q15303

Item Number: CM0015272



### Introduction

Rabbit polyclonal antibody targeting human ErbB4/HER4, validated for WB and ELISA, with cross-reactivity to rat. Ideal for research on neuregulin/EGF signaling, development, and oncogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | ErbB4/HER4 Rabbit pAb                                                                                        |
| Remarks/Alias    | HER4; ALS19; p180erbB4; ErbB4/HER4                                                                           |
| Species          | Human                                                                                                        |
| GeneID (Human)   | 2066                                                                                                         |
| GeneID           | 2066                                                                                                         |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1150-1250 of human ErbB4/HER4 (NP_005226.1) |
| Source           | Rabbit                                                                                                       |
| Category         | Polyclonal Antibodies                                                                                        |
| Application      | WB, ELISA                                                                                                    |
| Cross Reactivity | Rat                                                                                                          |
| SWISS            | Q15303                                                                                                       |
| Protein Weight   | 147kDa                                                                                                       |
| Shipping         | Ice bag                                                                                                      |