

Anti-Erc1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Iud2

Item Number: CM0025793

Introduction

Rabbit polyclonal antibody against human ERC1 (ELKS/Rab6IP2), validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying NF-kB signaling, synaptic organization, and membrane trafficking. Weight: 128 kDa.

[Learn More](#)

Parameter	Value
Product Name	ERC1 Rabbit pAb
Remarks/Alias	ELKS; Cast2; ERC-1; RAB6IP2; ERC1
Species	Human
GeneID (Human)	23085
GeneID	23085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-850 of human ERC1 (NP_829884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUD2
Protein Weight	128kDa
Shipping	Ice bag