

Anti-Erc2 Polyclonal Antibody For Wb, Ihc-P, Elisa - O15083

Item Number: CM0022300

Introduction

Rabbit polyclonal antibody against ERC2 (CAST/CAST1), validated for WB, IHC-P, ELISA with mouse and rat cross-reactivity. Ideal for studying presynaptic active zone organization and neurotransmitter release.

[Learn More](#)

Parameter	Value
Product Name	ERC2 Rabbit pAb
Remarks/Alias	CAST; CAST1; ELKSL; Spc110; SPBC110; ERC2
Species	Human
GeneID (Human)	26059
GeneID	26059
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 770-850 of human ERC2 (NP_056391.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15083
Protein Weight	111kDa
Shipping	Ice bag