

## Anti-Ercc1 Rabbit Monoclonal Antibody For Wb And Elisa - P07992

Item Number: CM0010424



### Introduction

Rabbit monoclonal antibody against human ERCC1 (P07992) validated for Western blot (WB) and ELISA. A key tool for studying DNA repair, interstrand crosslink repair, and genome stability.

[Learn More](#)

| Parameter        | Value                                                                                                        |
|------------------|--------------------------------------------------------------------------------------------------------------|
| Product Name     | ERCC1 Rabbit mAb                                                                                             |
| Remarks/Alias    | UV20; COFS4; RAD10; ERCC1                                                                                    |
| Species          | Human                                                                                                        |
| GeneID (human)   | 2067                                                                                                         |
| GeneID           | 2067                                                                                                         |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ERCC1 (P07992). |
| Source           | Rabbit                                                                                                       |
| Category         | Monoclonal Antibodies                                                                                        |
| Application      | WB, ELISA                                                                                                    |
| Cross Reactivity | Human                                                                                                        |
| SWISS            | P07992                                                                                                       |
| Protein Weight   | 33kDa                                                                                                        |
| Shipping         | Ice bag                                                                                                      |