

Anti-Ercc1 Rabbit Polyclonal Antibody For Wb And Elisa - P07992

Item Number: CM0026727



Introduction

Validate the anti-ERCC1 rabbit polyclonal antibody for WB and ELISA. Detects human, mouse, rat ERCC1, a key DNA repair protein. Ideal for research in genomic stability and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	ERCC1 Rabbit pAb
Remarks/Alias	UV20; COFS4; RAD10; ERCC1
Species	Human
GeneID (Human)	2067
GeneID	2067
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human ERCC1 (NP_973730.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07992
Protein Weight	33kDa
Shipping	Ice bag