


Anti-Ercc3 Rabbit Polyclonal Antibody For Wb, Elisa - P19447

Item Number: CM0015149



Introduction

High-quality rabbit polyclonal antibody against human ERCC3 (XPB), validated for WB and ELISA. Suitable for human, mouse, and rat samples. Ideal for transcription and DNA repair research.

[Learn More](#)

Parameter	Detail
Product Name	ERCC3 Rabbit pAb
Remarks/Alias	XPB; BTF2; Ssl2; TTD2; GTF2H; RAD25; TFIIH; ERCC3
Species	Human
GeneID (Human)	2071
GeneID	2071
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 512-782 of human ERCC3 (NP_000113.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19447
Protein Weight	89kDa
Shipping	Ice bag