

Anti-Ercc4 (Xpf) Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q92889

Item Number: CM0030085



Introduction

Rabbit polyclonal antibody against human ERCC4 (XPF) validated in WB, IP, ELISA; cross-reacts with mouse and rat. Targets DNA repair endonuclease XPF, essential for nucleotide excision repair, interstrand cross-link repair, and genome stability research.

[Learn More](#)

Parameter	Value
Product Name	ERCC4 Rabbit pAb
Remarks/Alias	XPF; RAD1; FANCD1; XFEPS; ERCC11; ERCC4
Species	Human
GeneID (Human)	2072
GeneID	2072
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 513-782 of human ERCC4 (NP_005227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92889
Protein Weight	104kDa
Shipping	Ice bag