


Anti-Ercc5 Polyclonal Antibody For If/Icc, Elisa - P28715

Item Number: CM0020513



Introduction

ERCC5 Rabbit Polyclonal Antibody validated for IF/ICC and ELISA. Detects human, mouse, rat ERCC5 (XPG), a key DNA repair endonuclease. Suitable for nucleotide excision repair and genomic stability research.

[Learn More](#)

Parameter	Value
Product Name	ERCC5 Rabbit pAb
Remarks/Alias	XPG; UVDR; XPGC; COFS3; ERCM2; ERCC5-201; ERCC5
Species	Human
GeneID (Human)	2073
GeneID	2073
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ERCC5 (NP_000114.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28715
Protein Weight	133kDa
Shipping	Ice bag