


Anti-Ercc6 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q03468

Item Number: CM0023563



Introduction

Rabbit polyclonal antibody against human ERCC6, validated for WB, IHC-P, ELISA with cross-reactivity in mouse and rat. Suitable for DNA repair and chromatin remodeling research.

[Learn More](#)

Parameter	Value
Product Name	ERCC6 Rabbit pAb
Remarks/Alias	CSB; CKN2; COFS; ARMD5; COFS1; POF11; RAD26; UVSS1; CSB-PGBD3; ERCC6
Species	Human
GeneID (Human)	2074
GeneID	2074
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 700-800 of human ERCC6 (NP_000115.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03468
Protein Weight	168kDa
Shipping	Ice bag