

Anti-Ercc6L2 Rabbit Polyclonal Antibody For Wb And Elisa - Q5T890

Item Number: CM0016165



Introduction

Rabbit polyclonal antibody targeting ERCC6L2 (RAD26L). Validated for WB and ELISA, cross-reacts with human and mouse. Aids DNA double-strand break repair and chromatin maintenance research. Predicted MW 177 kDa.

[Learn More](#)

Parameter	Value
Product Name	ERCC6L2 Rabbit pAb
Remarks/Alias	HEBO; BMFS2; SR278; RAD26L; C9orf102; ERCC6L2
Species	Human
GeneID (Human)	375748
GeneID	375748
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1262-1561 of human ERCC6L2 (NP_064592.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5T890
Protein Weight	177kDa
Shipping	Ice bag