

Anti-Ercc8 Rabbit Polyclonal Antibody For Wb, Elisa - Q13216

Item Number: CM0028887



Introduction

High-quality rabbit polyclonal antibody targeting human ERCC8 (CSA). Validated for Western blotting (WB) and ELISA. Cross-reacts with mouse. Ideal for research on transcription-coupled nucleotide excision repair, Cockayne syndrome, and DNA damage response. Predicted molecular weight: 44kDa.

[Learn More](#)

Parameter	Specification
Product Name	ERCC8 Rabbit pAb
Remarks/Alias	CSA, CKN1, UVSS2, ERCC8
Species	Human
GeneID (Human)	1161
GeneID	1161
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q13216
Protein Weight	44kDa
Shipping	Ice bag